

# RF Combo Jector

AUTO INJECTION & RF NEEDLING SYSTEM



The New Era of the  
**Cordless Auto Injector with RF**




#BIPOLAR  
#RF  
#INJECTOR

BIPOLAR

RADIO  
FREQUENCY

## Cordless Gun for Auto Injection with Bipolar RF during one shot



### Features



Medical CE Approved



RF · Injection · Microneedling **at once**




Compatible with **Fillers, Skin Boosters**





**Precise** Control : Depth, Dose, Duration

### 4 operating modes :

- 1) Injection + RF
- 2) Injection only
- 3) RF only
- 4) Filler ( 1-pin syringe )



## Preclinical test



Conducted by a dermatologist, Beom Jun Kim M.D.  
at Chung-Ang University Medical Center(Dept. of Dermatology)

## Acne Scar treatment (Age 23 / Male)

- RF Level : 4
- Frequency : 1
- Depth : 0.8 mm / 1.4 mm



Pre-Treatment

After 2 Months

After 5 Months

Conducted by dermatologist, In Kyu Lee M.D.

## RF Combo Jector **Principles**



**Treatment Area**

### Indications

- Acne Scars
- Pore
- Wrinkle
- Pigmentation
- Skin Rejuvenation


### Double Impact at once



STEP 01  
((o)) RF Treatment




STEP 02  
Suction & Injection




STEP 03  
Collagen Regeneration

### Various **needle cartridges**

#### • 01 RF & Injection needle



#### • 02 RF Treatment



#### • 03 Syringe




## Intuitive Interface & Customization

- Designed with basic buttons to prevent malfunction providing safe treatment environment
- Set the values per function, quick access available through short-cut key



## Powerful suction capabilities

- Pain relief by skin - adherence
- Minimizes loss of solution



## Specification & Composition

### Specification

<b>MODEL</b>	: RF Combo Jector ( Bipolar RF )
<b>DEPTH</b>	: 0.2mm ~ 5.0mm ( Max. )
<b>DIMENSION(mm)</b>	: W250 × D250 × H200 mm
<b>WEIGHT</b>	: 1.0kg
<b>INPUT</b>	: AC 100-240V, 50/60Hz, 0.3-0.5A
<b>BATTERY</b>	: Li-ion battery (3.7V / 3,000mAh)
<b>OPTIONAL</b>	: Auxiliary Cart

#### \* Auxiliary Cart

Cart is sold separately and not related to main unit`s functionality  
Designed for enhanced mobility and storage convenience

### Composition

