

rm

S T A R

★



Hybrid Skin Solution

Cordless RF & Microneedling Device



#RF

#MTS

#Needling

● ● ● Cordless RF & Microneedling Needling Device




Features

- **Medical CE Approved / FDA 510(K)TM** Cleared
- **Adjustable** Needle Length : 0 ~ 2 mm
- **Highly Effective in...**
Acne Scar Treatment & Wrinkle Reduction

1) The information provided applies to each model individually (RM STAR / RF STAR EX), as certification may vary between models.

2 operating Mode :

- RF & Microneedling at once
- Microneedling Only



Clinical test ● ● ●

Clinical Results



Wrinkle Treatment (Age 58 / Female)

RF Level : 3
Frequency : 3
Depth : 0.5 mm



Wrinkle Treatment (Age 62 / Female)

RF Level : 3
Frequency : 3
Depth : 0.3 mm



Acne Scar Treatment (Age 49 / Male)

RF Level : 5
Frequency : 3
Depth : 0.5 ~ 1.5 mm
(Around eyes / cheek)



Acne Scar & Pigmentation Treatment (Age 33 / Female)

RF Level : 3
Frequency : 3
Depth : 1.2 mm




Treatment Area

Indication


- Acne Scars
- Pore
- Wrinkle
- Pigmentation
- Skin Rejuvenation

Double Impact at once



Features


- Sterilized by EO Gas
- Medical Grade Material
- 24K gold plated
- Individual Blister Package



Handy but **powerful features**

Function

- Bipolar RF (1 Mhz)
- RF Level & Frequency : 1 ~ 5
- Depth : 0 ~ 2.0 mm
- Speed Level : 1 ~ 5
- Streamlined UI for easy and precise adjustment



Composition



Specification

MODEL	: RM STAR
ENERGY TYPE	: RF (Bipolar RF)
FREQUENCY	: 1 Mhz
DEPTH	: 0 mm ~ 2.0 mm
DIMENSION (mm)	: W 130 X D200 X H 125, Weight 0.5kg
INPUT	: AC 100 - 240V, 50 / 60Hz, 0.15 ~ 0.3A
BATTERY	: Li-ion battery, 700mAh, 3.1V



rm^{EX}
STAR
★

FDA 510(k) - Approved Model

- Equipped with an **additional LED main body** for a more **convenient treatment environment** for practitioners
- Approved for sales across **U.S.** and **globally**



* Cart is available for separate purchase



The Bom Building

322, GonghangDaero, Gangseo-gu, Seoul, Korea

T. +82 2. 523. 8295

F. +82 2. 523. 8297

E-mail

bomtech@bomtech.net

info@bomtech.net

global@bomtech.net

Website

www.bomtech.net

www.bomtechshop.com (Bomtech Shop)